




E-Xtreme

De esencia a extremo
Más pequeño | Más ligero | Más tranquilo
que nunca



Más tranquilo

Eighteeth motor innovador,
reduciendo el ruido desde el interior,
lo que permite experimentar la tranquilidad.



Más ligero



Menos de 100 gramos,
reduciendo la fatiga de la mano durante el tratamiento,
proporcionando la extrema experiencia de E-xtreme.

Más pequeño

- 360° rotación contra-ángulo



- Contra-ángulo en miniatura garantiza una mejor visibilidad y un acceso posterior más fácil



Contra-ángulo en miniatura
Super mini cabeza y cabello delgado



cabeza y cabello normal

- Diseño delgado y compacto para máxima comodidad y movilidad





E-xtreme

Make Essence to Extreme

Smaller

Lighter


Quieter than ever



Specifications

E-CAN Mini Contra Angle

Weight: 23.5 g





E-xtreme



99 g



120-650 rpm



0.5-4.0 N·cm



REC, FWD, REV



360°



800 mAh



E-connect Pro



160 g



120-1000 rpm



0.5-4.0 N·cm



REC, FWD, REV



340°



1500 mAh



Connect with E-PEX



E-value



160 g



120-1000 rpm



0.5-4.0 N·cm



ATC, REC, FWD, REV




360°



1900 mAh



E-connect S



170 g



120-1000 rpm



0.5-4.0 N·cm




ATC, REC, FWD, REV



340°



1500 mAh



Built-in Apex Locator